

# DENUDING micropipettes

Microtech Denuding micropipettes are prepared from borosilikate glass tubing (O.D. 1,2 mm, I.D. 0,75 mm), total length 9,5 cm have a blunt opening and are sold straight.

They are available in a choice of 3 sizes I.D. of opening: 300µm, 200µm and 150µm. The colour identity white (I.D. 120-150µm), grey (I.D. 180-200µm) and braun (I.D. 250-300µm).




4 pcs in packed

For single use only

40 pcs in box

## Identity colour

### WHITE



Article number

code

I.D.

005-120-C

005


120

005-150-C

005

150

### GREY



Article number

code

I.D.

005-180-B

005


180

005-200-B

005

200

### BROWN



Article number

code

I.D.

005-250-A

005

250

005-300-A

005

300